

Creating consistent comparative response curves of Cochlear Implant microphones



1. Cochlear Nucleus 6 CP910 BTE – DCTEST3 Lead only.

2. Cochlear Nucleus 7 CP1000 / Nucleus 8 CP1110 BTE – DCTEST4 lead, 3.5mm Mono plug to 3.5mm stereo socket adaptor, Monitor Adaptor Z544578.

3. Advanced Bionics Naida Q70 BTE – DCTEST4 lead and Naida Listening Check. CI-5823.

4. Advanced Bionics Naida CI M and Sky CI M – DCTEST4 lead and Advanced Bionics Marvel Listening Check CI-5827

5. Med-EL Rondo – DCTEST4 lead, Microphone Test Device 02883, Mini Battery Pack, Rondo Connection cable and MBPMTD C.Cable 08082.

6. Med-EL Sonnet BTE – DCTEST4 lead, Microphone Test Device 02883, Sonnet MTD Adaptor No.31545. Latest **version** of the FineTuner Remote. No.03132.

7. Med-EL Sonnet 2 BTE – DCTEST4 lead, Microphone Test Device 02883, Sonnet MTD Adaptor No.31545. Fine Tuner Echo Remote control No.37572*

* The FineTuner and Fine Tuner Echo are compatible with the Sonnets but may not work if user has paired the other remote. Check users remote control.